CARDIAC ASSIST

Manufactured by
Datascope Corp., 15 Law Drive - CN 40011,
Fairfield, New Jersey 07004-0011 U.S.A.
Tel: (973) 244-6100 Fax: (973) 244-6279

Authorized Representative:
Datascope B.V.
Drs. W. van Royenstraat 8
3871 AN Hoevelaken,
The Netherlands
Tel: +31.33.25.44.911 Fax: +31.33.25.44.393

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Datascope Medical Co. Ltd.
Lakeview Court
Ermine Business Park
Huntingdon
Cambridgeshire PE29 6XR
United Kingdom

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Datascope Corp.
Unit A, 30th Floor
Morrison Plaza
Morrison Hill Road
Wan Chai
Hong Kong

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CE 0044

CAUTION: Federal (U.S.A.) Law restricts this device to sale,
distribution and use by or on the order of a physician. Read this
entire booklet before opening inner pack


Fidelity is a trademark of Datascope Corp.
True Sheathless is a trademark of Datascope Corp.
Percor-STAT is a registered trademark of Datascope Corp.
Sheathless GARD is a trademark of Datascope Corp.
STAT-GARD is a registered trademark of Datascope Corp.

EXPLANATION OF SYMBOLS ON PACKAGE LABELING

	Do Not Reuse
REF	Catalog Reference Number
SN	Product Serial Number
	Lot/Batch Code
	Date of Manufacture (YYYY-MM-DD)
	Expiration Date (YYYY-MM-DD)
	Sterile. Ethylene Oxide Sterilization Method
	See Instruction for Use
	CE Conformity Marking in Accordance with the European Council Directive 93/42/EEC Concerning Medical Devices

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
PRODUCT STORAGE AND HANDLING REQUIREMENTS

Product sterile if packaging is undamaged or unopened.

DO NOT resterilize.

Handle with care. Do not expose to direct sunlight.

INTRA-AORTIC BALLOON (IAB) CATHETER

I. INDICATIONS FOR USE

- Refractory Unstable Angina.
- Impending Infarction.
- Acute Myocardial Infarction.
- Refractory Ventricular Failure.
- Complications of Acute MI (ie. Acute MR or VSD or papillary muscle rupture)
- Cardiogenic Shock.
- Support for diagnostic, percutaneous revascularization and interventional procedures.
- Ischemia related intractable ventricular arrhythmias.
- Septic Shock.
- Intraoperative pulsatile flow generation.
- Weaning from cardiopulmonary bypass.
- Cardiac support for non-cardiac surgery.
- Prophylactic support in preparation for cardiac surgery.
- Post-surgical myocardial dysfunction/low cardiac output syndrome.
- Myocardial Contusion.
- Mechanical bridge to other assist devices.
- Cardiac support following correction of anatomical defects.

II. CONTRAINDICATIONS

- Severe aortic insufficiency.
- Abdominal or aortic aneurysm.
- Severe calcific aorta-iliac disease or peripheral vascular disease.
- Introduction of the IAB catheter without the use of an introducer sheath is not recommended in patients with severe obesity, scarring of the groin or other contraindications to percutaneous insertion.

III. SUMMARY OF WARNINGS AND PRECAUTIONS

A WARNINGS

- If you continue to pump an IAB that has a leak, gaseous embolic injury of organs may result or a large blood clot may form within the balloon membrane requiring surgical removal of the IAB catheter.
- Do not inflate the IAB using a syringe or any other means if a balloon membrane leak is suspected.
- Perforation of a balloon membrane may indicate that the patient's vascular condition may induce abrasion or perforation in subsequent balloon membranes.
- The practitioner must be aware of adverse effects associated with percutaneous sheath introduction including bleeding at the insertion site, limb ischemia, infection, vessel trauma, and thrombosis.
- The practitioner must be aware of the potential for air embolism associated with open needles, sheaths, or catheter lumens in the patient's vasculature.
- Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body fluid precautions in the care of all patients.
- Do not cut the guide wire.
- Do not withdraw the guide wire against the needle bevel to avoid possibly severing or damaging the guide wire.

- If you encounter difficulty while inserting the IAB catheter without the use of an introducer sheath, remove the IAB catheter and insert the supplied introducer over the guide wire. From that point, continue with the remainder of the instructions to insert the IAB catheter with the use of an introducer sheath and guide wire.
- Use only the introducer dilator and sheath that are packaged with the insertion kit.
- Do not cut the introducer to alter the length.
- Blood leakage may occur if the hemostasis valve is removed from or not secured properly to the introducer sheath.
- Whenever possible, use fluoroscopy during IAB catheter insertion to ensure proper placement.
- Do not insert the IAB catheter unless the inner lumen is supported by a guide wire.
- Do not use excessive force when inserting the IAB catheter. If you use excessive force when inserting the IAB catheter, arterial tearing, dissection, or balloon membrane damage may occur.
- Any kinking or damage to the inner lumen may result in subsequent fatigue failure to the inner lumen when pumping.
- If you did not use fluoroscopy during IAB catheter insertion, you **MUST** take an X-ray as soon as possible to ascertain that the IAB catheter is properly positioned.
- If the balloon membrane has not entirely exited the introducer sheath, it will not inflate and deflate properly.
- Do not insert any portion of the Universal Sheath Seal below the skin line.
- NEVER INJECT AIR INTO THE INNER LUMEN** (female luer hub).
- Do not perform manual inflation of the IAB with the catheter extender tubing attached to the IAB catheter.
- If you aspirate blood from the male luer fitting of the extracorporeal tubing, immediate removal of the IAB catheter is indicated as damage may have occurred to the balloon membrane during insertion.
- If you note any restraint, unusual folding patterns of the balloon membrane or flexing of the inner lumen, immediately reposition the IAB catheter. Restraint may result in unusual folding patterns of the balloon membrane or flexing of the inner lumen for dual lumen IAB catheters. The life of a balloon membrane may be unpredictably shortened as a result of restraint which could lead to a leak.
- Do not attempt to withdraw the balloon membrane through the introducer sheath.
- If you feel any undue resistance during withdrawal of the IAB catheter, discontinue withdrawal and consider removal of the IAB catheter via an arteriotomy. Difficult removal may be a result of entrapment, due to a dried blood clot having formed within the balloon membrane from a balloon membrane leak.
- If limb ischemia is observed after IAB catheter removal, a vascular procedure may be indicated.
- Do not use excessive pressure when injecting contrast media through the inner lumen. Do not use an angiographic flow rate injector. The high pressure generated by the injector may damage the inner lumen. Do not use a syringe smaller than 20 cc to inject through the inner lumen. If resistance is met, consider the inner lumen obstructed and permanently seal it off.

B. Precautions

- Whenever possible, use fluoroscopy during guide wire and introducer sheath insertion.
- Use the appropriate sized guidewire for each Intra-aortic balloon catheter. Use only the .025" (.06 cm) guide wire with the 8 Fr. Fidelity IAB Catheter, the .030" (.08 cm) guide wire with the 9.5 Fr. True Sheathless IAB Catheter and the .030" (.08 cm) guide wire with the 10.5 Fr. Percor STAT-DL IAB Catheter.
- Use care to prevent kinking of the introducer during insertion.
- Pinching or kinking the reinforced introducer sheath may damage it, preventing insertion of the IAB catheter.
- For 8 Fr. IAB only, if the one-way valve is accidentally removed, cut the Gas Lumen Insert off at the base of the one-way valve to enable use of the one-way valve. Continue by following the instructions for use and disregard any further instructions regarding the Gas Lumen Insert.
- Do not remove IAB catheter from tray until immediately prior to inserting over guide wire.
- Remove the IAB catheter from the tray by pulling STRAIGHT out to avoid difficulty removing it from the tray or to avoid damaging it.
- Take care not to kink or place undue force on the IAB catheter.
- Always advance in short, continuous, one inch (2.5cm) strokes to avoid kinking the IAB Catheter.
- DO NOT twist the catheter during insertion.
- When using a chest X-ray to identify the position of the IAB Catheter, it is best to place the IABP in stand-by to improve visualization while the X-ray is taken and then immediately resume pumping.
- Do not place any sutures or ligatures around the outside diameter of the introducer to avoid kinking or damaging the sheath.
- For optimal signal quality, use no more than 8 feet (Maximum) of pressure tubing between the transducer and female luer hub of the Y-fitting.

- Use a standard flushing apparatus for arterial pressure monitoring with the inner lumen. Careful technique should be used in the set up and flushing of the arterial pressure monitoring apparatus to minimize the risk of an embolus entering the aorta where it could potentially enter the carotid or coronary arteries.
- Aspirate and discard a 3 cc volume of blood from the inner lumen prior to attaching a flushing apparatus to the female luer hub.
- Ensure that all air bubbles are removed from the inner lumen and flushing apparatus. In addition, tap the Y-fitting to remove all air bubbles.
- Prior to fast flushing, stop IAB pumping to reduce the risk of an embolus entering the aortic arch should an embolus be ejected from the inner lumen.
- For optimal signal quality, the inner lumen should not be used for blood sampling.
- Always aspirate 3 cc initially if the inner lumen aortic pressure line or the inner lumen becomes damped. If you meet resistance during aspiration, consider the inner lumen to be occluded. Discontinue the use of the inner lumen by placing a luer cap on the female luer hub.
- The use of in-line filters or other devices can potentially alter the appearance of the arterial pressure waveform.
- Do not over-tighten connections.
- The IAB should not remain inactive (i.e., not inflating and deflating) for more than 30 minutes because of the potential for thrombus formation.
- To avoid cutting the IAB catheter or the introducer, do not use scissors to remove the dressing.
- Inspect the entire introducer and IAB catheter upon removal to be certain the entire device has been removed.

IV. ADVERSE EFFECTS

A. Balloon Membrane Perforation

Balloon membrane perforation may be caused by:

- Contact with a sharp instrument.
- Fatigue failure due to unusual (biaxial) folding of the balloon membrane during use.
- Contact with calcified plaque resulting in abrasion of the surface and eventual perforation.

If perforation occurs, blood may be visible in the IAB catheter. If balloon membrane perforation is suspected, as may be evidenced by: 1) IAB pump leak alarms, 2) dried blood particles or serosanguineous fluid seen in the extracorporeal tubing or catheter extender or 3) a sudden change in the diastolic augmentation pressure waveform, the following procedure must be performed immediately:

- Stop pumping.
- Remove the IAB catheter
- Placing the patient in the Trendelenburg position should be considered if a leak is suspected.
- Consider IAB catheter replacement if the patient's condition warrants.

WARNINGS:

- If you continue to pump an IAB with a leak, gaseous embolic injury of organs may result or a large blood clot may form within the balloon membrane requiring surgical removal of the IAB catheter.
- Do not inflate the IAB using a syringe or any other means if a balloon membrane leak is suspected.
- Perforation of a balloon membrane may indicate that the patient's vascular condition may induce abrasion or perforation in subsequent balloon membranes.

NOTE: The length of time a balloon membrane can survive such contact with plaque or unusual folding is unpredictable. A leak in an IAB catheter within the bloodstream may allow gas to enter the patient's bloodstream which may result in patient injury. Large perforations are rare, therefore the small quantity of gas released is usually asymptomatic. The rate of incidence at each individual hospital may be influenced by the degree of vascular disease in that patient population, by the location of the IAB catheter in the aorta, or by using a balloon membrane size inappropriate for the specific patient.

B. Limb Ischemia

During or after IABP, limb ischemia may result. It can be caused by an obstruction of flow due to:

- Thrombus formation.
- Creation of an intimal layer separation or flap.
- The presence of the introducer sheath or IAB catheter.

After IAB catheter removal, if limb ischemia is observed, a vascular procedure may be indicated. Monitor distal limb for the development of compartment syndrome.

C. Bleeding at the Insertion Site

Bleeding at the insertion site may be caused by:

- Trauma to the artery during insertion of the IAB.
- Excessive catheter movement at the insertion site.
- Anticoagulation.

Bleeding at the insertion site may be controlled with direct pressure at the insertion site, assuring adequate distal blood flow. If bleeding persists, surgical repair of the insertion site may be indicated.

D. Infection

Infection may occur due to interruption of normal skin integrity at the IAB catheter insertion site. Sterile technique should be used during IAB catheter insertion and during dressing changes.

Assess the patient for the development of IAB catheter related infection and treat, if necessary.

E. Thrombocytopenia

Thrombocytopenia may develop due to mechanical damage to the platelets. Monitor platelet count and replace platelets, if necessary.

F. Aortic Dissection

Aortic dissection may occur during insertion of the IAB catheter. Symptoms can include back and/or abdominal pain, a decreased hematocrit, and hemodynamic instability.

G. Thrombosis

Thrombus formation may occur during counterpulsation. The symptoms associated with thrombosis formation and treatment will depend on the organ system involved.

V. EQUIPMENT REQUIRED

The following is a list of compatible consoles for the Datascope IAB catheters at heart rates not exceeding 140 BPM. The use of Datascope catheters on some consoles at higher heart rates may result in decreased augmentation. Consoles at higher heart rates may result in decreased augmentation.

IAB Catheters	Datascope Systems	Arrow Pumps
Fidelity 8 Fr. 9.5 Fr. True Sheathless 10.5 Fr. Percor STAT-DL	98, 98XT and CS100	Arrow ACAT, AutoCat

IAB catheter sizing

Select the most appropriate size IAB catheter for the patient from the following balloon membrane sizing table:

Balloon Membrane Volume (cc)	Balloon Membrane Dimensions		Patient Height	
	Length (mm)	Diameter (mm)	(ft)	(cm)
25	174	14.7	<5'	<152
34	219	14.7	5'0"-5'4"	152-162
40	263	15	5'4"-6'0"	162-183
50	269	16.3	>6'0"	>183

NOTE: This information is to be used only as a guideline. Clinical judgment and patient factors (i.e., torso length) should be considered when selecting the most appropriate size IAB catheter.

The following sterile equipment is required for insertion of the IAB catheter. Inspect all components prior to use.

The following items are supplied by Datascope and are provided with all products:

- One sterile one-way valve and 60 cc syringe (included with IAB catheter)
- One sterile Insertion Kit which includes:
 - One 18 gauge angiographic needle
 - One vessel dilator
 - One introducer sheath with hemostasis valve
 - One introducer dilator
 - One three-way stopcock
 - One luer cap
 - One Datascope Catheter Extender
 - Pressure Tubing:
 - One 4 ft. (1.2 m) length of pressure tubing (with 9.5 Fr. and 10.5 Fr. IABs)
 - Two 4 ft. (1.2 m) lengths of pressure tubing (with 8 Fr. IABs)
 - Guide Wires:
 - 8 Fr. IAB – One .025" (.06 cm) x 145 cm PTFE coated guide wire
 - 9.5 Fr. IAB – One .030" (.08 cm) x 145 cm PTFE coated guide wire
 - 10.5 Fr. IAB – One .030" (.08 cm) x 145 cm PTFE coated guide wire

Optional item: One sterile Arrow Catheter Extender

The following items are not supplied by Datascope:

- Local anesthetic with syringe and needle
- One sterile scalpel and blade
- Radiographic contrast material
- One sterile 20 cc syringe
- Sterile lint free sponges

VI. INSTRUCTIONS

WARNINGS:

- The practitioner must be aware of adverse effects associated with percutaneous sheath introduction including bleeding at the insertion site, limb ischemia, vessel trauma, and thrombosis.
- The practitioner must be aware of the potential for air embolism associated with open needles, sheaths, or catheter lumens in the patient's vasculature.
- Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.

PRECAUTION: Whenever possible, use fluoroscopy during guide wire and introducer sheath insertion.

A. IAB CATHETER INSERTION

NOTE: Insertion technique is for Sheathless Insertion of the IAB unless otherwise indicated. The Percor STAT-DL 10.5 Fr. IAB must be inserted using the Introducer Sheath.

- Make the customary preparations for percutaneous catheterization and administer appropriate local anesthesia.
- Insert the angiographic needle into the common femoral artery at a 45 degree or less angle. (See Figure 1)


Figure 1

- Insert the J-tip end of the guide wire through the angiographic needle and advance into the thoracic aorta.

WARNING: Do not cut the guide wire.

WARNING: Do not withdraw the guide wire against the needle bevel to avoid possibly severing or damaging the guide wire.

PRECAUTION: Use the appropriate sized guidewire for each Intra-aortic balloon catheter. Use only the .025" (.06 cm) guide wire with the 8 Fr. Fidelity IAB Catheter, the .030" (.08 cm) guide wire with the 9.5 Fr. True Sheathless IAB Catheter and the .030" (.08 cm) guide wire with the 10.5 Fr. Percor STAT-DL IAB Catheter.

- Keeping the guide wire in place, remove and discard the needle.
- Wipe the blood from the guide wire with a wet, lint-free sponge.
- Make a small incision at the exit of the guide wire to facilitate inserting the vessel dilator through the skin.
- Place the tapered end of the vessel dilator over the exposed guide wire and dilate the artery by pushing the vessel dilator into the arterial lumen.
- Keeping the guide wire in place, remove and discard the vessel dilator. Apply pressure at the wound site to control bleeding.
- Wipe the blood from the guide wire with a wet, lint-free sponge.
- For sheathless insertion: Spread the tissue at the incision with a tissue dilator.

WARNING: If you encounter difficulty while inserting the IAB catheter using sheathless insertion, remove the IAB catheter and insert the supplied introducer sheath over the guide wire. From that point, continue with the remainder of the instructions to insert the IAB catheter with the use of an introducer sheath.

A1. INSTRUCTIONS FOR INSERTING THE INTRODUCER SHEATH

- Make the customary preparations for percutaneous catheterization and administer appropriate local anesthesia.
- Insert the angiographic needle into the common femoral artery at a 45 degree or less angle. (See Figure 1)
- Insert the J-tip end of the guide wire through the angiographic needle and advance into the thoracic aorta.
- Keeping the guide wire in place, remove and discard the needle.

WARNING:

- Do not cut the guide wire.
- Do not withdraw the guide wire against the needle bevel to avoid possibly severing or damaging the guide wire.
- Use only the introducer dilator and sheath that are packaged with the insertion kit.
- Do not cut the introducer to alter the length.
- When using the 8Fr. sheath, insert the introducer dilator into the sheath hub and twist to lock in place. When using the 10 Fr. or 11.5 Fr. introducer, assure that the dilator is secured to the sheath with the locking nut.

WARNING: Blood leakage may occur if the hemostasis valve is removed from or not secured properly to the introducer sheath.

- Place the tapered end of the introducer dilator over the exposed guide wire and advance the introducer with a rotary motion into the arterial lumen. (See Figure 2)
- Remove the dilator leaving the sheath in place.
- Continue from Step 11 of Section A: IAB CATHETER INSERTION.


Figure 2

PRECAUTIONS:

- Use care to prevent kinking of the introducer during insertion.
- Pinching or kinking the reinforced introducer sheath may damage it, preventing insertion of the IAB catheter.


NOTE: During insertion of the IAB catheter, after removal of the introducer dilator, some blood leakage may be observed past the hemostasis valve. This blood leakage will subside as the IAB catheter is advanced into the introducer sheath. Do not kink or pinch the reinforced introducer sheath to control bleeding.

WARNINGS:

- Whenever possible, use fluoroscopy during IAB catheter insertion to ensure proper placement.
- Do not insert the IAB catheter unless the inner lumen is supported by a guide wire.
- It is recommended that anticoagulant therapy be administered in accordance with standard hospital practice.

NOTE: Maintain vacuum on the IAB catheter throughout insertion. Do not remove the one-way valve.

- Remove the tray from the sterile packaging.
- Firmly attach the one-way valve to the male luer fitting of the extracorporeal tubing. (See Figure 3)

PRECAUTION:

If the one-way valve is accidentally removed, cut the Gas Lumen Insert off at the base of the one-way valve to enable use of the one-way valve. Continue by following the instructions for use and disregard any further instructions regarding the Gas Lumen Insert.


figure 3

- With the 60 cc syringe, slowly aspirate at least 30 cc (See Figure 3). Remove the syringe while leaving the one-way valve in place.
- Remove the extracorporeal tubing portion of the IAB from the tray – (DO NOT disconnect one-way valve when removing the extracorporeal tubing from the tray.)


figure 4

PRECAUTIONS:

- Do not remove IAB catheter from tray until immediately prior to inserting over guide wire.
 - Remove the IAB catheter from the tray by pulling STRAIGHT out to avoid difficulty removing it from the tray or to avoid damaging it.
 - Take care not to kink or place undue force on the IAB catheter.
- Remove the IAB catheter by lifting the Y-fitting and catheter from the tray.
 - Withdraw the balloon membrane from the protective sleeves by PULLING THE IAB CATHETER STRAIGHT OUT OF THE SLEEVES THROUGH THE SLOT PROVIDED IN THE TRAY. (See Figure 4.) **THE SLEEVES AND RETAINING LABEL ARE NOT PART OF THE IAB CATHETER AND SHOULD REMAIN ATTACHED TO THE TRAY.**
 - If fluoroscopy is not used, measure the distance from the angle of Louis or between the second and third intercostal space extending down to the umbilicus and then obliquely over to the femoral insertion site. Slide the sheath seal or Sheathless GARD protector up the catheter to this measured distance.
 - Lubricate the IAB catheter by passing it through a basin of sterile saline solution.
 - Remove the stylet wire from the inner lumen (Figure 5). Do not attempt stylet reinsertion.
 - Manually flush inner lumen with 3-5 cc flush solution.


figure 5

- Insert guide wire through the inner lumen (See Figure 6.) Advance the IAB catheter over the guide wire until the guide wire exits the female luer hub. Always ensure the operator has complete control of the guidewire.


figure 6

WARNING: Do not use excessive force when inserting the IAB catheter. If you use excessive force when inserting the IAB catheter, arterial tearing, dissection, or balloon membrane damage may occur.

NOTE: During insertion of the IAB catheter, arterial blood under pressure may run down the length of the folds in the balloon membrane and drip or be expelled under arterial pressure from the balloon membrane/catheter junction. THIS "CHANNELING" IS NOT A LEAK. As the IAB catheter is advanced the bleeding will diminish.

- While controlling the proximal end of the guide wire, advance the IAB catheter over the guide wire into the artery. **Always grasp the IAB catheter no more than one inch (2.5cm) from the insertion site or sheath hub and advance in short continuous strokes to avoid kinking the IAB catheter while maintaining complete control of the guidewire.**

PRECAUTIONS:

- Always advance in short, continuous, one inch (2.5cm) strokes to avoid kinking the IAB Catheter.
- DO NOT twist the catheter during insertion.

WARNINGS:

- Any kinking or damage to the inner lumen may result in subsequent fatigue failure to the inner lumen when pumping.
 - Do not insert the IAB catheter unless the inner lumen is supported by a guidewire.
- Advance the IAB catheter to the proper position in the descending thoracic aorta, with the IAB catheter tip just distal (approximately 2cm) to the left sub-clavian artery. (See Figure 7)


Figure 7

WARNING: If you do not use fluoroscopy during IAB catheter Insertion, you MUST take an X-ray as soon as possible to ascertain that the IAB catheter is properly positioned.

Precaution: When using a chest X-ray to identify the location of the IAB Catheter, it is best to place the IABP in stand-by to improve visualization while the X-ray is taken and then immediately resume pumping.

A2. WHEN USING INTRODUCER SHEATH

WARNING: If the balloon membrane has not entirely exited the introducer sheath, it will not inflate and deflate properly.

NOTE: Either during advancing or once the catheter is correctly positioned; confirm that the IAB catheter membrane has fully exited the sheath. See Figure 8.




figure 8

- The first single band from the IAB catheter tip indicates that the entire balloon membrane has exited from the 6" (15.2cm) introducer sheath/hemostasis valve and may now be inflated.
- Single bands follow at 3/4" (1.9cm) increments.
- A double band indicates that the entire balloon membrane has exited from the 11" (27.9cm) introducer sheath/hemostasis valve (10.5 Fr. IAB's only) and may now be inflated.
- Maintain sterility of the exposed catheter until the proper IAB catheter position is verified.
- After proper placement of the IAB Catheter, push the Universal Sheath Seal as close to the insertion site as possible. (See Figures below)

NOTE: If blood is seen passing the sheath seal following insertion through a sheath, disengage sheath seal from hemostasis valve.



WARNING: Do not insert any portion of the Universal Sheath Seal below the skin line.


WARNING: Do not insert the Sheathless GARD protector into the subcutaneous tissue.

A3. INSTRUCTIONS FOR USE OF SHEATH SEAL WITH INTRODUCER SHEATH

PRECAUTION: Do not place any sutures or ligatures around the outside diameter of the Introducer to avoid kinking or damaging the sheath.

Push the sheath seal down over the hemostasis valve hub.

Push the sheath seal down over the hemostasis valve hub.

If it becomes necessary to reposition the IAB catheter, hold the sheath seal in one hand, grasp the catheter through the STAT-GARD sleeve with the other hand, and reposition under aseptic conditions. Do not attempt to reposition the IAB catheter by moving the sheath.

B. INITIATING IAB PUMPING (IABP)

NOTE: Do not raise the head of the bed greater than 45°.

- After positioning the IAB catheter, remove the guide wire.
- Once the catheter is in place, aspirate and discard 3 cc of blood from the inner lumen and then immediately perform a manual flush using a syringe filled with 3 cc to 5 cc of flush solution. This will minimize the chances of stagnant blood clotting in the inner lumen.

NOTE: If you are not using the inner lumen for pressure monitoring, aspirate 3 cc of blood from the inner lumen and cap the female luer hub with the supplied luer cap.

WARNING: NEVER INJECT AIR INTO THE INNER LUMEN (female luer hub).

- Release the vacuum from the IAB catheter by removing the one-way valve from the IAB catheter male luer.
- For Fidelity 8 Fr. IAB's only, after IAB insertion has been completed, remove the Gas Lumen Insert by loosening the one-way valve from the extracorporeal tubing's male luer and slowly withdraw the Gas Lumen Insert from the gas passage of the IAB catheter (See Figure 9).


Figure 9

- Connection of the IAB catheter to the pump.

NOTES:

- Ensure that all connections are leak-free.
- All catheter extenders are sterile and should only be used one time.
- Use one catheter extender when connecting the IAB catheter to the IAB pump.

- When a Datascope IAB pump is used:** Connect the IAB catheter's male luer fitting to the female luer fitting of the catheter extender. Connect the male luer fitting of the catheter extender to the safety chamber/disk.
- When an Arrow IAB pump is used:** Connect the IAB catheter's male luer fitting to the female luer fitting of the Arrow catheter extender, then connect the catheter extender to the appropriate system. Adjust the volume setting on the Arrow pump, according to the operating instructions, to match the IAB catheter volume.

NOTE: Prior to using the Arrow catheter extender, remove the protective tubing from the pins.

- Follow the IAB pump operating instructions to initiate pumping. If augmentation is not within the desired range, see Appendix A, Factors Affecting Augmentation.
- If after a few cycles of counterpulsation, it appears that the balloon membrane is not fully open, perform the following procedure:

WARNING: Do not perform manual inflation of the IAB with the catheter extender tubing attached to the IAB catheter.

- Detach the catheter extender from the IAB catheter's male luer fitting.
- Attach the supplied 3-way stopcock and syringe to the IAB catheter's male luer fitting.

WARNING: NEVER INJECT AIR INTO THE INNER LUMEN (female luer hub).

- Aspirate to assure blood is not returned through the extracorporeal tubing.

WARNING: If you aspirate blood from the male luer fitting of the extracorporeal tubing, immediately remove the IAB catheter as damage may have occurred to the balloon membrane during insertion.

- Inflate the IAB with air or helium, as follows, and IMMEDIATELY aspirate:

IAB Catheter Volume	Inflate Volume
25 cc.	50 cc.
34 cc., 40 cc or 50 cc.	60 cc.

- Remove the 3-way stopcock and syringe and reattach the IAB catheter's male luer fitting to the catheter extender and resume pumping.
- Ensure that the balloon membrane is inflating and deflating in an unrestrained manner and not restrained as a result of lodging beneath a plaque, within a subintimal space, within the subclavian artery, aortic arch, abdominal aorta, or if the balloon membrane volume is too large for the specific patient's aorta.

WARNING: If you note any restraint, unusual folding patterns of the balloon membrane or flexing of the inner lumen, immediately reposition the IAB catheter. The life of a balloon membrane may be unpredictably shortened as a result of restraint, which could lead to a balloon failure.

- If you observe unusual bleeding or subcutaneous hematoma at the insertion site, treat appropriately.
- Evaluate peripheral pulses. If the distal pulse is not satisfactory or signs of limb ischemia are present, exercise discretion concerning the continuation of IAB pumping.
- Restrain movement of the IAB catheter by securing the suture pads and the Y-fitting to the skin using a StatLock Securement Device (see Figure 10) or sutures. You may also secure the Y-fitting to the patient with tape. (See Figure 11)
- Apply a dressing to the insertion site using sterile technique according to hospital policy.


figure10


figure11

C. PRESSURE MONITORING THROUGH THE INNER LUMEN


Figure 12

NOTE: ALL FITTINGS MUST BE LEAK FREE

PRECAUTION: For optimal signal quality, use no more than 8 feet (Maximum) of pressure tubing between the transducer and female luer hub of the Y-fitting.

When monitoring pressure through the inner lumen, use a standard arterial pressure monitoring apparatus connected to a three-way stopcock. (See Figure 12.) Connect the three-way stopcock to the female luer hub of the inner lumen. A 3 cc/hour continuous flow through the inner lumen is recommended. The anticoagulation dosage should be in accordance with standard hospital practice for arterial pressure lines and may be modified, on physician discretion, for patients receiving anticoagulation therapy. Per hospital policy, a fast forward flush may be performed hourly to help maintain patency of the inner lumen.

PRECAUTIONS DURING PRESSURE MONITORING THROUGH IAB CATHETER

- Use a standard flushing apparatus for arterial pressure monitoring with the inner lumen. Careful technique should be used in the set up and flushing of the arterial pressure monitoring apparatus to minimize the risk of an embolus entering the aorta where it could potentially enter the carotid or coronary arteries.
- Aspirate and discard a 3 cc volume of blood from the inner lumen prior to attaching a flushing apparatus to the female luer hub.
- Ensure that all air bubbles are removed from the inner lumen and flushing apparatus. In addition, tap the Y-fitting to remove all air bubbles.
- Prior to fast flushing, stop IAB pumping to reduce the risk of an embolus entering the aortic arch should an embolus be ejected from the inner lumen.
- For optimal signal quality the inner lumen should not be used for blood sampling.
- Always aspirate 3 cc initially if the inner lumen aortic pressure line or the inner lumen becomes damped. If you meet resistance during aspiration, consider the inner lumen to be occluded. Discontinue the use of the inner lumen by placing a luer cap on the female luer hub.
- The use of in-line filters or other devices can potentially alter the appearance of the arterial pressure waveform.

- Do not over-tighten connections.

RECOMMENDATIONS FOR ACHIEVING OPTIMAL PRESSURE SIGNAL QUALITY

- Use no more than 8 ft. (2.4 m) of a low compliance pressure tubing such as that supplied by Datascope in the IAB Insertion Kit between the transducer and Y-fitting of the catheter.
- Once the catheter is in place, aspirate and discard 3 cc of blood from the inner lumen and then immediately perform a manual flush using a syringe filled with 3 cc to 5 cc of flush solution. This will minimize the chances of stagnant blood clotting in the inner lumen.
- Apply only gentle force to the syringe when aspirating the inner lumen.
- Do not use a R.O.S.E. (Resonance Over Shoot Eliminator) or other damping device.
- Remove air from flush bag prior to pressurizing.
- Prime the pressure set-up using gravity flush.
- Maintain 300 mmHg of pressure on the flush solution and elevate it above the transducer.
- Whenever the inner lumen of the IAB becomes filled with blood (such as after aspiration), the flush valve should be activated for a **minimum of 15 seconds** in addition to the time it takes to clear the pressure tubing of blood.
- Ensure that all air bubbles are removed from the inner lumen and flushing apparatus.
- Use room temperature flush solution.

D. MANUALLY INFLATING AND DEFLATING THE IAB CATHETER

PRECAUTION: The IAB catheter should not remain inactive (i.e., not inflating and deflating) for more than 30 minutes because of the potential for thrombus formation.

To keep the IAB catheter active in the event of pump failure, manually inflate and deflate the IAB as follows:

WARNING: Do not perform manual inflation of the IAB with the catheter extender tubing attached to the IAB catheter.

- Detach the catheter extender from the IAB catheter's male luer fitting.
- Attach the supplied 3-way stopcock and syringe to the IAB catheter's male luer fitting.

WARNING: NEVER INJECT AIR INTO THE INNER LUMEN (female luer hub).

- Aspirate to assure blood is not returned through the extracorporeal tubing.

WARNING: If you aspirate blood from the male luer fitting of the extracorporeal tubing, immediate removal of the IAB catheter is indicated as damage may have occurred to the balloon membrane during insertion.

- Inflate the IAB with 40 cc air or helium and immediately aspirate. Repeat every 5 minutes while the IAB is inactive.
- Remove the 3-way stopcock and syringe and reattach the IAB catheter's male luer fitting to the catheter extender and resume pumping.

E. REMOVING THE IAB CATHETER

- Consider tapering or discontinuing anticoagulation therapy prior to removal.
- Stop IAB pumping.
- Disconnect the IAB catheter from the IAB pump permitting the IAB catheter to vent to atmosphere. Patient blood pressure will collapse the balloon membrane for withdrawal.
- Remove all sutures and dressings.

PRECAUTION: To avoid cutting the IAB catheter or the Introducer, do not use scissors to remove the dressing.

- Remove the IAB Catheter.

- If an introducer sheath is used: Loosen the sheath seal from the hub and withdraw the IAB catheter through the introducer sheath until resistance is met.
- Remove the IAB catheter and the introducer sheath as a unit

WARNING: If you feel any undue resistance during withdrawal of the IAB catheter, discontinue withdrawal and consider removal of the IAB catheter via an arteriotomy. Difficult removal may be a result of entrapment, due to a dried blood clot having formed within the balloon membrane from a balloon membrane leak.

- Apply digital pressure below the puncture site during IAB catheter removal. Allow free proximal bleeding for a few seconds, then apply pressure above the puncture site and allow a few seconds of back bleeding. Establish hemostasis to the puncture.
- Carefully examine the limb distal to the insertion site for adequate perfusion.

WARNING: If limb ischemia is observed after IAB catheter removal, a vascular procedure may be indicated.

PRECAUTION: Inspect the entire Introducer sheath and IAB catheter to be certain the entire device has been removed.

NOTE: If you must initiate IAB pumping after the removal of the IAB catheter, you can perform a percutaneous insertion on the contralateral femoral artery. Do not use the same insertion site.

APPENDIX A: FACTORS AFFECTING AUGMENTATION

If after pumping has commenced, augmentation is not within the desired range, one of the following may be indicated:

- Balloon membrane has not fully exited from the introducer sheath.** Pull the introducer sheath back until the balloon membrane fully exits the introducer sheath.
- Balloon membrane has not fully opened.** See Section B, Initiating IAB Pumping (IABP), Item 7.
- IAB augmentation/volume control on the IAB pump is set too low.** Adjust the IAB augmentation/volume control on the IAB pump.
- IAB catheter is positioned in the aortic arch, subclavian artery, or otherwise malpositioned in the aorta.** Observe the IAB catheter under fluoroscopy. If malpositioned, remove any suture that may have been applied over the sheath seal or STAT-GARD sleeve and reposition the IAB catheter.
- IAB catheter is positioned in a false lumen.** Use fluoroscopy to check for the appropriate intra-luminal position of the IAB catheter tip, as follows:
 - Aspirate a 3 cc volume of blood from the inner lumen of the IAB catheter (female luer hub).
 - With the IAB pumping, inject 10-20 cc of contrast media through the inner lumen of the IAB catheter.

WARNING: Do not use excessive pressure when injecting contrast media through the inner lumen. Do not use an angiographic flow rate injector. The high pressure generated by the injector may damage the inner lumen. Do not use a syringe smaller than 20 cc to inject through the inner lumen. If resistance is met, consider the inner lumen obstructed and permanently seal it off.

- If the IAB catheter tip is appropriately positioned within the aorta, the radiopaque cloud will wash away with the next two or three heart beats.
- If the IAB catheter tip is positioned in a false lumen, the radiopaque cloud will remain surrounding the IAB catheter.
- If the IAB catheter is found to lie within a false lumen, remove the IAB catheter from the patient. Consider insertion of a new IAB catheter into the contralateral femoral artery.

In addition to the above, physiological conditions can contribute to poor augmentation. Among them are:

- Patient's mean arterial blood pressure is low.
- Patient's systemic vascular resistance is low.
- Patient's heart rate is rapid enough to compromise ventricular filling and ejection.

APPENDIX B: LIMITED WARRANTY

Datascope Corp. warrants that all its IAB catheters are free from defects in workmanship and materials for a period of three years from the date of purchase or until the expiration date is reached or until the date of use, whichever comes first. Datascope Corp. shall not be liable for any incidental, special or consequential loss, damage, or expense directly arising from the use of this product. Liability under this warranty and the buyer's exclusive remedy is limited to replacement of the product which, under normal use and service, shall have been found by the Company to be defective in materials or workmanship. It shall be the buyer's obligation to return any such product to the Company for examination for replacement liability.

No agent, employee, or representative of Datascope Corp. has any authority to bind Datascope Corp. to any affirmation, representation, or warranty concerning its products, and any affirmation, representation, or warranty made by any agent, employee or representative shall not be enforceable by the buyer.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ANY OTHER EXPRESS OR IMPLIED WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS, AND OF ANY OTHER OBLIGATION ON THE PART OF THE SELLER.

Damage to any product or parts through misuse, neglect, accident, or by affixing any nonstandard accessory attachments or by any customer modification voids this warranty. Datascope Corp. makes no warranty whatever in regard to trade accessories, such being subject to the warranty of their respective manufacturers.

A condition of this warranty is that this equipment or any accessories which are claimed to be defective be returned, when authorized by Datascope, freight prepaid, to Datascope Corp., 15 Law Drive, CN 40011, Fairfield, New Jersey 07004-0011. Datascope Corp. shall not have any responsibility in the event of loss or damage in transit.

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